

透過低頻率的全身律動誘導神經突處的Pc12m3細胞生長因子增生以修復受損的神經突觸

▶▶▶ Low-Frequency, Whole Body Vibration Induced Neurite Outgrowth by Pc12m3 Cells with Impaired Nerve Growth Factor-Induced Neurite Outgrowth

發表自：

Journal of Novel Physiotherapies, 2015, 5:1
新型物理治療學期刊, 2015年02月12日, 第05卷, 第01期

簡述內容：

隨著年齡、壓力、疾病、感染、失調等狀態可能會出現神經失衡或反應不敏銳現象，當此現象出現亦可能造成後續例如顫抖、疼痛、無知覺、肌肉萎縮等不同情況，尤其像嚴重疾病如中風、脊椎損傷患者這類無法自行活動者更加速神經退化，而老年者的神經敏感度不佳對後續生活亦造成嚴重危害風險。

本研究證實以低頻率(7Hz)的全身垂直律動，即可誘使神經細胞增生、刺激神經纖維生長、促進神經連結，遠比單獨神經生長因子(NGF)的誘導高於3倍，表示全身垂直律動是促使神經增生的有效方式，後續運用於神經問題方面、衰弱、疾病或無法自主運動者，將是最佳改善的訓練方式之一。



Low-Frequency, Whole Body Vibration Induced Neurite Outgrowth by PC12m3 Cells with Impaired Nerve Growth Factor-Induced Neurite Outgrowth

Yoshihisa Koike^{1*}, Reiko Tutida¹, Yuko Hayashi¹, Yoko Yamanishi¹ and Yoshio Kano²

¹Department of Occupational Therapy, Faculty of Health and Welfare, Prefectural University of Hiroshima, Mihara City, Hiroshima 723-005, Japan

²Department of Occupational Therapy, School of Health Science, Kibi International University, 8 Iga-machi Takahashi City, Okayama 716-8508, Japan

*Corresponding author: Yoshihisa Koike, Department of Occupational Therapy, Faculty of Health and Welfare, Prefectural University of Hiroshima, Mihara City, Hiroshima 723-005, Japan, Tel: +81-848-60-1212; E-mail: koike@pu-hiroshima.ac.jp

Received: October 30, 2014, Accepted: February 05, 2015, Published: February 12, 2015

Copyright: © 2015 Koike Y, et al. This is an open-access article distributed under the terms of the Creative Commons Attribution License, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited.

Abstract

Purpose: To investigate the effects of very low-frequency whole body vibration (WBV) stimulation using PC12 m3 cells.

Methods: Mutant, drug-hypersensitive, PC12 m3 cells were obtained by means of continuous culture of neural PC12 cells and then stimulated by low-frequency vibration at frequencies of 7–13 Hz for 5–20 min with simultaneous exposure to nerve growth factor (NGF).

Results: The frequency of neurite outgrowth by PC12 m3 cells induced by 7-Hz vibratory stimulation was approximately 3-fold greater than that induced by NGF alone. Moreover, activated cyclic-AMP responsive element (CRE)-binding protein (CREB) expression was induced in PC12 m3 cells stimulated by 7-Hz vibration.

Conclusion: Low-frequency vibratory stimulation induces neurite outgrowth via a p38 mitogen-activated protein kinase/CREB signaling pathway in PC12 m3 cells. Together, these results suggest that WBV is an efficient method to stimulate neurite growth.

Keywords: Whole body vibration; p38 MAPK; CREB; PC12m3 cell

Introduction

Recent studies have reported that whole body vibration (WBV) was effective for improving bone strength, muscle strength in the lower extremities, and functional mobility [1-3]. However, the vibratory frequency ranges that may be beneficial have not been standardized. In particular, the effects of very low frequency vibrations are uncertain. Thus, the purpose of our study was to determine the effects of very low-frequency vibratory stimulation to induce neurite outgrowth by PC12m3 cells.

Our basic investigation focused on a non-pharmacological therapy for the elderly, which resulted in fewer adverse events as compared with pharmacological therapy regimens [4]. The effects of each treatment were examined using PC12 m3 cells before the initiation of the clinical study [5,6]. The treatments found to be successful at this level were further studied clinically. In the clinical study, Vibroacoustic therapy was found to provide relaxation effects for elderly nursing home residents, and it improved symptoms of depression [7]. The "steam foot spa" treatment resulted in improved cognitive function and a temporary decrease in high blood pressure in geriatric inpatients [8]. The effect of this treatment was confirmed in PC12 m3 cells and further observed in ongoing clinical investigations.

For this purpose, we used a device that could generate very low-frequency vibrations (Tap Master). Vibrations with this device could be artificially produced to generate vertical sinusoidal vibrations. Moreover, the amplitudes were in a low frequency range of 5.6–13 Hz.

With this device, we investigated the effects of WBV stimulation on neurite outgrowth by a drug-hypersensitive PC12 mutant cell line, PC12m3 cells. PC12 cells are para-neuronal cells of rat adrenal pheochromocytoma origin. PC12 cells are a useful undifferentiated cell model that is sensitive to the activity of nerve growth factor (NGF), which induces their differentiation to nerve cells and neurite extensions [9,10]. Sustained activation of mitogen-activated protein kinase (MAPK) plays an important role in neurite outgrowth by PC12 cells [11]. In addition, because of the widespread use of PC12 cells under various culture conditions, spontaneous variants often arise [12,13]. We developed a drug-hypersensitive cell line, PC12m3 cells, for which neurite outgrowth could be induced by various drugs, such as FK506, cAMP, and Calcimycin [14].

Mammalian cells have at least three MAPK pathways that regulate the activities of extracellular signal-regulated kinase (ERK; also known as p42 and p44 MAPK), c-Jun N-terminal kinase (JNK), and p38 mitogen-activated protein kinase (p38 MAPK). The ERK pathway is important for regulating cellular differentiation. JNK pathways have been implicated in regulating neuron apoptosis. However, the p38 MAPK signaling pathway's specific system has not been elucidated [15]. In this study, we found that activated p38 MAPK was involved in the mechanism by which neurite outgrowth by PC12m3 cells was induced by very low-frequency vibration. We also found that low-frequency vibration induced the expression of activated cyclic-AMP response element (CRE)-binding protein (CREB). CREB is a transcription factor that is a target for MAPK.

Materials and Methods

Reagents and cell culture

Nerve growth factor (NGF; 2.5S) was from Takara (Osaka, Japan). The p38 MAPK inhibitor SB203580 was from Sigma (St. Louis, MO). PC12m3 cells were cultured in Dulbecco's modified Eagle's medium (DMEM) supplemented with 0.35% glucose, 10% horse serum, 5% fetal bovine serum (FBS), and 100 units/ml of kanamycin. Cells were grown at 37°C in a 5% CO₂ atmosphere for all experiments described below.

Whole body vibration device

A whole body vibration device (Tap Master, Bodygreen Health Co., Ltd Taiwan) was used to generate both sinusoidal wave vibrations and vertical wave vibrations. The frequency range of this device was 5.6 Hz to 13 Hz and the amplitude range was -1 mm to 1 mm. Our aim in these experiments was to identify differences in the effects of vibration at 3-Hz intervals as compared with those at 7-, 10-, and 13-Hz intervals.

Determination of neurite outgrowth

A single-cell suspension of PC12m3 cells was obtained by trituration in DMEM. For cell culture experiments, cells were plated at a density of 2-5×10⁵ cells per plate on 25 cm² plates with serum-containing DMEM, after which cells were immediately treated with NGF and different very low-frequency vibrations. After incubation for 5 to 7 days, the frequency of neuritogenesis was determined by measuring neurite lengths and numbers. Cells that had one or more neurites with lengths >1.5-fold of the cell body diameter were measured as previously described [15,16]. Results were expressed as the means ± S.E.M.'s for 100-200 cells sampled in three independent experiments.

Analysis for activated p38 MAPK and CREB induction

Induction of activated p38 MAPK and CREB expression was determined as described previously [17]. In brief, PC12m3 cells were plated at a density of 1×10⁶ cells/25 cm² in a flask with serum-containing medium and cultured for 3 days. Next, culture medium was replaced with 0.5% FBS-containing medium for 48 h. Prior to NGF treatment, cells were incubated for 2 h in serum-free medium. Cells were then stimulated for 5 to 20 min by adding NGF (30 ng/ml) only or by vibratory stimulation at 7 Hz to 13 Hz along with NGF (30 ng/ml). Activated p38 MAPK and CREB expression was then assessed using cell lysates.

Cells were lysed in lysis buffer. Lysate aliquots (10-15 µg) for each sample were separated on an SDS 10% polyacrylamide gel and then transferred to polyvinylidene difluoride membranes. Membranes were probed with antibodies directed against phospho-p38 MAPK or phospho-CREB (New England BioLabs; Beverly, MA) used at a dilution of 1:1000 in blocking buffer (5% nonfat dry milk) at 4°C for 12 h. The blots were then reacted with a secondary antibody, horseradish peroxidase-conjugated anti-rabbit IgG (diluted 1:2000 in blocking buffer) at room temperature for 60 min. Blots were stained for 1 min using a chemiluminescence reagent (LumiGLO chemiluminescent reagent, Kirkegaard and Perry Laboratories) and exposed to x-ray film.

Data analysis

Results were analyzed by one-way analysis of variance (ANOVA) and Dunnett's test to identify significant differences from the control. The significance level was chosen as 0.05.

Results

Establishing PC12m3 cells

During continuous culture of neural PC12 cells, we obtained a drug-hypersensitive PC12 mutant cell line, PC12m3 cells. PC12 cells are subject to spontaneous mutations that result in generating PC12 variants. The variant cell line that we obtained had impaired NGF-induced neurite outgrowth. When these cells were cultured for 2 weeks under acidic conditions of Cl⁻ several surviving clones appeared. By using the ring isolation procedure, ten colonies (designated PC12m, PC12m, and so on) were selected and propagated in mass culture.

PC12m3 cells subsequently exhibited poor neurite outgrowth in response to NGF. PC12m3 cells that were treated with NGF showed enhanced neurite outgrowth in response to various agents, such as FK506, c-AMP, and Calcimycin. Furthermore, these cells exhibited sustained activation of ERK that was induced by different agents [15].

Vibration induced neurite outgrowth by PC12m3 cells

PC12m3 cells were tested for their sensitivity to physical stimulation using vibratory stimuli. Figure 1 shows photomicrographs of PC12m3 cells after treatment with NGF alone and those that were treated with low-frequency WBV stimuli at 7 Hz, 10 Hz, or 13 Hz for 10 min and in the presence of NGF. Vibratory stimulation at 7 Hz combined with NGF resulted in enhanced neurite outgrowth as compared to NGF treatment alone. However, vibratory stimulation at 13 Hz had no effect on neuritogenesis by PC12m3 cells (Figure1).




Figure 1: PC12m3 cell neurite outgrowth induction by low-frequency vibratory stimulation. PC12m3 cells were exposed to vibratory stimulation immediately after NGF treatment (A, B, C). (A) PC12m3 cells treated with NGF alone, (B) PC12m3 cells stimulated with 7-Hz vibratory stimulation for 10 min and with NGF treatment, and (C) PC12m3 cells stimulated with 13-Hz vibratory stimulation for 10 min and with NGF treatment. Phase contrast photomicrographs of PC12m3 cells were acquired at 7 days after treatment (×200).

PC12m3 cells were then exposed to low-frequency WBV of different frequencies ranging from 7 to 13 Hz for 5 to 20 min along with NGF treatment (Figure 2). This showed that low-frequency vibration at 7 Hz induced enhanced neurite outgrowth, whereas

vibratory stimuli of 10 Hz and 13 Hz had no effect on neurite outgrowth. The frequency of neurite outgrowth induced by 7-Hz low-frequency vibration was approximately 3-fold greater than that induced by NGF alone.




Figure 2: Frequencies of neurite outgrowth by PC12m3 cells induced by vibratory stimulation ranging from 5.6 to 13 Hz for 5 to 20 min. Results are means \pm S.E.Ms (n = 3). By post-hoc Bonferroni correction analysis, the frequency of neurite outgrowth induced by 7-Hz low-frequency vibratory stimulation was significantly different from that of controls. *P<0.05.

Induction of activated p38 MAPK and CREB expression in PC12m3 cells by 7-Hz low-frequency vibration

Because activation of p38 MAPK plays an important role in neuron differentiation by PC12 cells [10], we examined whether low-frequency vibratory stimulation induced neurite outgrowth by PC12m3 cells was due to its effects on p38 MAPK activity. PC12m3 cells were stimulated with vibrations of 7 Hz and 13 Hz for 10 min or were not stimulated (controls), after which activated p38 MAPK (phospho-p38) expression was assessed by immunoblotting (Figure 3).




Figure 3: Induction of activated p38 MAPK and CREB expression in PC12m3 cells by 7-Hz low-frequency vibratory stimulation. PC12m3 cells were serum-starved and then stimulated or not by vibration for 30 min along with NGF treatment. After treatment, cells were lysed, and protein extracts were analyzed by Western blot using anti-phospho-p38 MAPK and anti-phospho-CREB (Ser 133) antibodies.

This showed that 7 Hz vibration had enhanced the expression of activated p38 MAPK. In addition, neurite outgrowth was inhibited by a specific p38 MAPK inhibitor, SB203580 (Figure 4). In contrast, activated p38 MAPK expression was not enhanced in cells stimulated

with 13 Hz vibration and in control cells. Thus, it appeared that low-frequency vibration induced neurite outgrowth involved the p38 MAPK signaling pathway in PC12m3 cells.

We also examined vibratory stimulation-induced CREB activation in PC12m3 cells. Cells were stimulated with vibration of 7 Hz, after which activated CREB (phospho-CREB) expression was assessed by immunoblotting (Figure 3). This showed that activated CREB expression was enhanced by vibratory stimulation of 7 Hz for 10 min.




Figure 4: Inhibition of neurite outgrowth by PC12m3 cells induced by low-frequency vibratory stimulation and with NGF treatment. PC12m3 cells were treated with NGF alone, stimulated with 7-Hz low-frequency vibratory stimulation for 10 min and with NGF treatment, or stimulated with 7-Hz low-frequency vibratory stimulation for 10 min with NGF and treated with a specific p38 MAPK inhibitor, SB203580 (2 μ l, 5 μ l). The percentages of cells with neurite outgrowth were determined after 7 days. *P<0.05.**P<0.01.

Discussion

Our results showed that VAT provided relaxation effects to elderly nursing home residents and promoted passive aerobic exercise, consequently improving symptoms of depression and improving sleep patterns [7].

The results of previous VAT studies using PC12 m3 cells exposed to vibratory stimuli at frequencies of 10–200 Hz for 30 min followed by NGF treatment showed that vibrations of 20–100 Hz (low frequency) enhanced neurite outgrowth. In the present study, the frequency of neurite outgrowth induced by 40-Hz (low frequency) vibrations was approximately three-fold greater than that induced by NGF alone. Activation of the p38 MAPK (mitogen-activated protein kinase) pathway plays an important role in neuronal differentiation of PC12 m3 cells. Therefore, we examined whether the ability of a low-frequency vibratory stimulus to induce neurite outgrowth of PC12 m3 cells is a reflection of its effect on p38 MAPK activity [5]. The results showed that vibratory stimuli induced neurite outgrowth via the p38 MAPK signaling pathway in PC12 m3 cells.

In studies on intracellular transduction pathways of aerobic exercise using mouse muscle cells, Akimoto proposed a pathway where MKK3/MKK6 in MAPKK is activated by aerobic exercise, leading to downstream p38MAPK activation and eventual biosynthesis of mitochondria [18]. Moreover, Falempin et al. suggested that tendon vibration (120 Hz) of murine soleus muscle can be used as a paradigm to counteract the atrophic process observed after hind-limb unloading

[19]. In addition, Blumenthal reported that aerobic exercise is effective for alleviating depression in the elderly, with an effect comparable to that of antidepressant drugs [20]. In consideration of these findings, mitigation of depression due to VAT for elderly NH residents seems to be enhanced by vibro-tactile stimuli due to passive aerobic exercise. Still, the p38MAPK signaling pathway has been reported to have protective effects on cardiac muscle [21].

Vascular endothelial cell growth factor (VEGF) induces endothelial cell proliferation and movement, remodeling of the extracellular matrix, the formation of capillary tubules, and vascular leakage [22]. Moreover, VEGF plays a crucial role in the development of the cardiovascular system and in promoting angiogenesis that is associated with physiological and pathological processes [23]. VEGF is an endothelial cell-specific mitogen that promotes numerous other phenomena that are necessary for angiogenesis [24].

Mitogen-activated protein kinase (MAPK) is important in the induction of endothelial cell proliferation that is induced by VEGF [25]. VEGF promotes the phosphorylation and activation of the cyclic AMP-responsive element binding protein (CREB) by activating p38 MAPK/mitogen- and stress-activated protein kinase-1 (MSK-1) signaling pathways that are downstream of the kinase insert domain receptor (KDR; vascular endothelial cell growth factor receptor-2) [26,27]. Because CREB is a transcription factor that plays an important role in promoting cellular proliferation and adaptive responses, a mechanism has been defined by which VEGF/KDR may allow endothelial cells to respond to changes in their environment through alterations in their gene expression [24]. These findings indicate that by mediating cellular responses to VEGF, CREB, through its activation of p38 MAPK/ MSK-1 signaling pathways, is likely to play an important role in endothelial cell function and blood development.

In the present study, p38 MAPK was activated by vibratory stimulation at 7 Hz for 10 min and promoted neurite outgrowth by PC12m3 cells of approximately 3-fold greater than that induced by NGF alone. Furthermore, neurite outgrowth induction by low-frequency vibration was inhibited by a specific p38 MAPK inhibitor, SB203580, in the presence of NGF. A recent study demonstrated that p38 MAPK participated in preserving neurite growth cones during neurite outgrowth and in regulating cellular differentiation and survival. These results suggest that p38 MAPK was not involved in neurite outgrowth that was induced by 7 Hz vibratory stimulation.

Moreover, we detected an activated cyclic-AMP responsive element (CRE)-binding protein (CREB) expression in PC12m3 cells after vibratory stimulation at 7 Hz. CREB is a transcription factor that is the target of a various signaling pathways that mediate cell responses to extracellular stimuli, including proliferation, differentiation, and adaptive responses. p38 MAPK can promote CREB binding to CRE. Our results suggested that p38 MAPK may induce neurite outgrowth by PC12m3 cells via a CREB signaling pathway. This indicates that p38MAPK may induce neurite outgrowth via a CREB signaling pathway after low-frequency whole body vibration.

Another recent study suggested that mechanical stimulation by whole body vibration increased shear stress at the walls of blood vessels, which resulted in increased blood flow velocity after vibration was terminated and promoted muscle deoxygenation [28,29]. Increased shear stress in capillaries in skeletal muscle can initiate capillary growth by producing a disturbance in the luminal side of the basement membrane [30]. It appears likely that shear stress can induce angiogenesis [31].

All of these findings indicate that whole body vibration induced CREB expression by activating a p38 MAPK signaling pathway. This signaling pathway is likely to play an important role in endothelial cell function and blood development.

Prospective Studies

These results suggest that WBV can effectively improve depression and heart disease, to an extent, in the elderly.

Acknowledgement

This work was support by Project for the Promotion of Medicine-Technological Industries collaboration in Hiroshima Prefecture

References

1. Lai CL, Tseng SY, Chen CN, Liao WC, Wang CH, et al. (2013) Effect of 6 months of whole body vibration on lumbar spine bone density in postmenopausal women: a randomized controlled trial. *Clin Interv Aging* 8: 1603-1609.
2. Marin PJ, Rhea MR (2010) Effects of vibration training on muscle strength: a meta-analysis. *J Strength Cond Res* 24: 548-556.
3. Delecluse C, Roelants M, Verschueren S (2003) Strength increase after whole-body vibration compared with resistance training. *Med Sci Sports Exerc* 35: 1033-1041.
4. Koike Y, Yamanishi Y, Kano Y (2014) Effects of Nonpharmacological Therapies for Diseases of the Elderly. *Psychology Research* 4: 389-396.
5. Koike Y, Iwamoto S, Kimata Y, Nohno T, Hiragami F, et al. (2004) Low-frequency vibratory sound induces neurite outgrowth in PC12m3 cells in which nerve growth factor-induced neurite outgrowth is impaired. *Tiss.Cult.Res. Commun* 2: 81-90.
6. Kano Y, Nakagiri S, Nohno T, Hiragami F, Kawamura K, et al. (2004) Heat shock induces neurite outgrowth in PC12m3 cells via the p38 mitogen-activated protein kinase pathway. *Brain Res* 1026: 302-306.
7. Koike Y, Hoshitani M, Tabata Y, Seki K, Nishimura R, et al. (2012). Effects of vibroacoustic therapy on elderly nursing home residents with depression. *J. Phys. Ther. Sci.* 24: 291-294.
8. Koike Y, Kondo H, Kondo S, Takagi M, Kano Y (2013) Effect of a steam foot spa on geriatric inpatients with cognitive impairment: a pilot study. *Clin Interv Aging* 8: 543-548.
9. Green LA (1977) A quantitative bioassay for nerve growth factor (NGF) activity employing a clonal pheochromocytoma cell line. *Brain Res* 133: 350-353.
10. Greene LA, Tischler AS (1976) Establishment of a noradrenergic clonal line of rat adrenal pheochromocytoma cells which respond to nerve growth factor. *Proc Natl Acad Sci USA* 73: 2424-2428.
11. Erhardt P, Troppmair J, Rapp UR, Cooper GM (1995) Differential regulation of Raf-1 and B-Raf and Ras-dependent activation of mitogen-activated protein kinase by cyclic AMP in PC12 cells. *Mol Cell Biol* 15: 5524-5530.
12. Hatanaka H (1981) Nerve growth factor-mediated stimulation of tyrosine hydroxylase activity in a clonal rat pheochromocytoma cell line. *Brain Res* 222: 225-233.
13. Teng KK, Georgieff IS, Aletta JM, Nunez J, Shelanski ML, et al. (1993) Characterization of a PC12 cell sub-clone (PC12-C41) with enhanced neurite outgrowth capacity: implications for a modulatory role of high molecular weight tau in neurogenesis. *J Cell Sci* 106: 611-626.
14. Kano Y, Hiragami F, Kawamura K, Kimata Y, Nakagiri S, et al. (2002) Immunosuppressant FK506 induces sustained activation of MAP kinase and promotes neurite outgrowth in PC12 mutant cells incapable of differentiating. *Cell Struct Funct* 27: 393-398.
15. Morooka T, Nishida E (1998) Requirement of p38 mitogen-activated protein kinase for neuronal differentiation in PC12 cells. *J Biol Chem* 273: 24285-24288.

16. Yao R, Osada H (1997) Induction of neurite outgrowth in PC12 cells by gamma-lactam-related compounds via Ras-MAP kinase signaling pathway independent mechanism. *Exp Cell Res* 234: 233-239.
17. Inoue S, Motoda H, Koike Y, Kawamura K, Hiragami F, et al. (2008) Microwave irradiation induces neurite outgrowth in PC12m3 cells via the p38 mitogen-activated protein kinase pathway. *Neurosci Lett* 432: 35-39.
18. Akimoto T, Pohnert SC, Li P, Zhang M, Gumbs C, et al. (2005) Exercise stimulates Pgc-1alpha transcription in skeletal muscle through activation of the p38 MAPK pathway. *J Biol Chem* 280: 19587-19593.
19. Falempin M, In-Albon SF (1999) Influence of brief daily tendon vibration on rat soleus muscle in non-weight-bearing situation. *J Appl Physiol* (1985) 87: 3-9.
20. Blumenthal JA, Babyak MA, Moore KA, Craighead WE, Herman S, et al. (1999) Effects of exercise training on older patients with major depression. *Arch Intern Med* 159: 2349-2356.
21. Nishida K, Yamaguchi O, Hirotani S, Hikoso S, Higuchi Y, et al. (2004) p38alpha mitogen-activated protein kinase plays a critical role in cardiomyocyte survival but not in cardiac hypertrophic growth in response to pressure overload. *Mol Cell Biol* 24: 10611-10620.
22. Ferrara N (1999) Molecular and biological properties of vascular endothelial growth factor. *J Mol Med (Berl)* 77: 527-543.
23. Soker S, Kaefer M, Johnson M, Klagsbrun M, Atala A, et al. (2001) Vascular endothelial growth factor-mediated autocrine stimulation of prostate tumor cells coincides with progression to a malignant phenotype. *Am J Pathol* 159: 651-659.
24. Mayo LD, Kessler KM, Pincheira R, Warren RS, Donner DB (2001) Vascular endothelial cell growth factor activates CRE-binding protein by signaling through the KDR receptor tyrosine kinase. *J Biol Chem* 276: 25184-25189.
25. Wu LW, Mayo LD, Dunbar JD, Kessler KM, Baerwald MR, et al. (2000) Utilization of distinct signaling pathways by receptors for vascular endothelial cell growth factor and other mitogens in the induction of endothelial cell proliferation. *J Biol Chem* 275: 5096-5103.
26. Deak M, Clifton AD, Lucocq LM, Alessi DR (1998) Mitogen- and stress-activated protein kinase-1 (MSK1) is directly activated by MAPK and SAPK2/p38, and may mediate activation of CREB. *EMBO J* 17: 4426-4441.
27. Tan Y, Rouse J, Zhang A, Cariati S, Cohen P, et al. (1996) FGF and stress regulate CREB and ATF-1 via a pathway involving p38 MAP kinase and MAPKAP kinase-2. *EMBO J* 15: 4629-4642.
28. Beijer A, Rosenberger A, Bölck B, Suhr F, Rittweger J, et al. (2013) Whole-body vibrations do not elevate the angiogenic stimulus when applied during resistance exercise. *PLoS One* 8: e80143.
29. Kerschlan-Schindl K, Grampp S, Henk C, Resch H, Preisinger E, et al. (2001) Whole-body vibration exercise leads to alterations in muscle blood volume. *Clin Physiol* 21: 377-382.
30. Yamada E, Kusaka T, Miyamoto K, Tanaka S, Morita S, et al. (2005) Vastus lateralis oxygenation and blood volume measured by near-infrared spectroscopy during whole body vibration. *Clin Physiol Funct Imaging* 25: 203-208.
31. Hudlicka O (1991) What makes blood vessels grow? *J Physiol* 444: 1-24.